



# *surgiTrend*

SURGEON'S CHAIR



ENJOY QUALITY AND DESIGN

# surgiTrend



Item No.: 2200 0012 with Comfort seat (Illustration)  
 Item No.: 2200 0010 with Dynamic seat

For back protection and maximum work flexibility in confined spaces. Sophisticated, space-saving drive mechanism with reduced size for ample foot space combines tilt stability with easy and smooth maneuvering. Outstanding sitting comfort and a wide range of adjusting functions allow every surgeon to find his or her personal best and healthiest working position for any work scenario. Available in two seating system types: Comfort or Dynamic.



## Seating systems



**Comfort**

Seat angle adjustment



Back angle adjustment



Back height adjustment



### Upholstery / colour card

Standard colour

schwarz  
 Item No.: 3037 0223



"Periline"

Page 111

**Comfort**



**Dynamic**

Dynamic-function



Lumbar support angle adjustment




Lumbar support height adjustment



Standard colour

schwarz  
 Item No.: 3037 0021



"Stamskin Top"

Page 110

**Dynamic**

## Armrests

- Stable, ergonomically shaped armrests with ample support surface (350 x 120 mm)
- Height-adjustable and equipped with a ball joint for 360-degree rotation
- Designed for personalized adjustment to the surgeon
- Optionally available with linear guide and one-hand release



## Chassis

- 4 high-quality, smooth-running double castors
- 5th castors for added resistance to tipping
- Optimized foot space
- Central braking system effective on 4 castors
- Integrated foot keypad panel for seat-height adjustment
- Power supply-independent battery box with emergency stop button



## Technical specifications

### Lifting column

- Drive 3,000 N (IPX4)
- Lifting height 200 mm

### Chassis

- W x L (incl. braking system/foot keypad) 600 x 650 mm
- Castors ø 100 mm
- 5th castors ø 65 mm
- Central brake acting on all four castors
- Emergency stop button

### Battery unit

- Battery box with 2 batteries (lead-gel), pollution free, maintenance-free, cadmium-free 2 x 12 V, 7.2 Ah
- Universal battery charger: Input: 100 – 240 V  
Output: 24 V

### Seating system: Comfort (Dynamic)

- Min. sitting height 550 mm (600 mm)
- Max. sitting height 750 mm (800 mm)
- Seat angle adjustments +5° / -5° (-12°)
- Back angle adjustments +5° / -25° (-16°)

### Weight

45 kg

### Max. load

150 kg

Subject to technical modifications.

## Technical dimensions

### Seating system Comfort



### Seating system Dynamic



# Components



*surgiTrend*  
Dynamic seat Item No.: 2200 0010



*surgiTrend*  
Comfort seat Item No.: 2200 0012



*surgiLine*  
Formed seat Item No.: 2200 0008



*surgiLine*  
Comfort seat Item No.: 2200 0013



Backup battery

Item No.: 3047 0000



*surgiTrend*  
Footswitch plate

Item No.: 3351 1011



*surgiLine*  
Footswitch plate

Item No.: 3351 0014



Armrest pair with  
linear guide and  
one-hand release,  
factory-fitted

Item No.: 3351 006



Armrest pair with  
one-hand release,  
factory-fitted

Item No.: 3351 1004



# "Stamskin Top"



User friendly: High flexibility of the material ensures a perfect finish.

Safe for the user: Phthalate-free, easy to clean and fulfills fire safety standards.

## Product Specification

- Extremely durable
- Phthalate free
- Easy to clean
- Flame protection
- Mildew resistant

*surgi*Trend  
*assist*Trend

Dynamic seat 2200 0010  
2500 0002

+ Item No.: 3037 0...



F4340-10120  
Item No.: .....023



F4340-07422  
Item No.: .....016



F4340-20130  
Item No.: .....084




F4340-20156  
Item No.: .....207



F4340-50622  
Item No.: .....215




F4340-50623  
Item No.: .....216




F4340-20129  
Item No.: .....146




F4340-00002  
Item No.: .....021



F4340-20162  
Item No.: .....205



F4340-20161  
Item No.: .....206



F4340-20150  
Item No.: .....204



F4340-5015  
Item No.: .....198



F4340-07496  
Item No.: .....046



F4340-07432  
Item No.: .....002



F4340-20154  
Item No.: .....208



F4340-07445  
Item No.: .....030



F4340-50621  
Item No.: .....199



F4340-07480  
Item No.: .....039



F4340-07436  
Item No.: .....020



F4340-10295  
Item No.: .....044



F4340-07459  
Item No.: .....028



F4340-20157  
Item No.: .....211



F4340-5347  
Item No.: .....209



F4340-50454  
Item No.: .....210



F4340-20149  
Item No.: .....212



F4340-50619  
Item No.: .....214



F4340-20127  
Item No.: .....079



F4340-07477  
Item No.: .....036



F4340-07476  
Item No.: .....035



F4340-50618  
Item No.: .....220



F4340-07485  
Item No.: .....034



F4340-50620  
Item No.: .....213



F4340-07478  
Item No.: .....037



F4340-07479  
Item No.: .....038

# "Periline"



## Product Specification

- Without Phthalates
- Biocompatibility:  
Cytotoxicity (EN ISO 10 993-5) and  
test for irritation and delayed-type  
hypersensitivity (EN ISO 10 993-10)
- Urine resistance
- Blood resistance
- Resistance to perspiration
- Resistance to spitting
- Light fastness
- Color fastness to rubbing
- Flame retardant

<i>surgiForce</i>	Comfort seat	2600 0004
<i>surgiForce light</i>	Comfort seat	2600 0006
<i>surgiTrend</i>	Comfort seat	2200 0012
<i>surgiLine</i>	Comfort seat	2200 0013
<i>assistPro F</i>		2400 0005
<i>assistPro H</i>		2400 0009
<i>easyPro F</i>		2400 0010
<i>easyPro H</i>		2400 0006
<i>easyPro R</i>		2400 0013

+ Item No.: 3037 0...



kastanie 015817  
Item No.: .....240



light brown 016037  
Item No.: .....239



creme 015812  
Item No.: .....225



weiß 010881  
Item No.: .....224



stein 015813  
Item No.: .....236



saphir 015814  
Item No.: .....231



delft 010766  
Item No.: .....241



antistatisch schwarz  
Item No.: .....223



schwarz 011499  
Item No.: .....235



anthracit 016046  
Item No.: .....234



azur 015815  
Item No.: .....230



aubergine 016041  
Item No.: .....229



lavendel 016042  
Item No.: .....228



rot 010878  
Item No.: .....238



orange 016040  
Item No.: .....237



kiwi 015816  
Item No.: .....233



lime green 016043  
Item No.: .....232



lemon 016039  
Item No.: .....226



safran 010879  
Item No.: .....227

The colours displayed on your device may differ slightly from the original colours. If you wish to receive original colour charts, please contact us at [info@ufsk-osys.com](mailto:info@ufsk-osys.com).

# Upgrade your OR to perfection.



*eyeForce*



*500 XLE*



*500 XLE comfort*



*500 bariatric*



*500 XLS*



*500 ECO*



*400 XL*



*motionLine*



*motionLine-E*



*surgiForce*



*surgiForce light*



*surgiTrend*



*surgiLine*



*stayClean wipes*



*stayClean*



*staySoft*



*stayClean roll*



*assistTrend*



*assistPro F*



*assistPro H*



*easyPro F*



*easyPro H*



*easyPro R*



*assistLine*



*supportLine*



*supportLine plus*

## CONTACT US



**UFSK**  
INTERNATIONAL

**OSYS**

UFSK-International OSYS GmbH  
Kirchhoffstraße 1  
93055 Regensburg  
Germany

Phone +49 (0)941 78862 15  
Fax +49 (0)941 78862 35  
Mail [info@ufsk-osys.com](mailto:info@ufsk-osys.com)  
Web [ufsk-osys.com](http://ufsk-osys.com)

**MADE IN GERMANY**  
Design | Construction | Manufacturing